




## Product Specification

Page 1

**Suporte para tubos redondos para Tubos redondos e podem ser combinados com os suportes de parede euroMatic® (95.963.001), Suporte de parede twin plus (95.963.007) e Suporte de parede vario 2000 (95.963.020)**



### Product description

Order number	95.963.006
Descrição do artigo	Suporte para tubos redondos, para Tubos redondos e podem ser combinados com os suportes de parede euroMatic® (95.963.001), Suporte de parede twin plus (95.963.007) e Suporte de parede vario 2000 (95.963.020), material: POM, (CxLxA): 130 x 40 x 88 mm

### Product characteristics

Indicado para	Tubos redondos e podem ser combinados com os suportes de parede euroMatic® (95.963.001), Suporte de parede twin plus (95.963.007) e Suporte de parede vario 2000 (95.963.020)
---------------	---

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

Page 2

**Suporte para tubos redondos para Tubos redondos e podem ser combinados com os suportes de parede euroMatic® (95.963.001), Suporte de parede twin plus (95.963.007) e Suporte de parede vario 2000 (95.963.020)**

### Size

Largura do produto	40 mm
Altura do produto	88 mm
Comprimento do produto	130 mm

### Material & colours

Material do produto	Polioximetileno (POM)
Cor do produto	prata/branca

### Purity & certification

Categoria do produto	Não é um dispositivo médico   DIV
Em lotes	não

### Packaging

Quantidade mínima do pedido	1
Tipo de embalagem interior	caixa de papelão
Unidades por caixa interna	1
Unidades por caixa	10
Unidades por palete	1680
Profundidade da embalagem externa	288 mm
Largura da embalagem externa	218 mm
Altura da embalagem externa	150 mm
Volume da embalagem externa	0.0094 cbm
Peso do produto	0.4003 kg
Peso da embalagem externa	4.23 kg
EAN da caixa interna	4038917364300
EAN da embalagem externa	4038917364294

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**