




Product Specification

**Microvette® CB 300 Heparina de lítio LH, 300 µl,
tampa laranja, tampa de pressão, fundo plano**

Page 1



Product description

| | |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number | 16.434 |
| Descrição do artigo | Microvette® CB 300 Heparina de lítio LH, Medicina veterinária, preparo: Heparina de lítio, volume nominal: 300 µl, (CxØ) com tampa: 46,5 x 10,8 mm, tampa laranja, tampa de pressão, fundo plano, recipiente interno: em forma de trompete, tubo: marrom, 100 unid./pacote |

Product characteristics

| | |
|------------------------------|------------------------|
| Tipo de retirada | capilar |
| Tipo de preparo | Heparina de lítio |
| Código de cor | EU |
| Tampa | Tampa macia com argola |
| Tipo de fechamento | tampa de pressão |
| Campo de aplicação | Medicina veterinária |
| Forma do fundo | fundo plano |
| Centrifugação recomendada | 2000 U/min |
| Tempo de centrifugação | 5 Min |
| Forma do recipiente interior | em forma de trompete |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Microvette® CB 300 Heparina de lítio LH, 300 µl,
tampa laranja, tampa de pressão, fundo plano**

Page 2

Size

| | |
|-----------------------|---------|
| Volume nominal | 300 µl |
| Diâmetro | 10.8 mm |
| Comprimento com tampa | 46.5 mm |
| Comprimento sem tampa | 43.7 mm |

Material & colours

| | |
|------------------------|----------------------------------|
| Material do produto | Polipropileno (PP) |
| Cor do produto | transparente |
| Material de fechamento | Low Density Polyethylene (LD-PE) |
| Cor da tampa | laranja |

Purity & certification

| | |
|----------------------|-----------------------------------|
| Categoria do produto | Não é um dispositivo médico DIV |
| Em lotes | sim |

Packaging

| | |
|--------------------------------------|---------------|
| Quantidade mínima do pedido | 1000 |
| Tipo de embalagem interior | pacote |
| Unidades por caixa interna | 100 |
| Unidades por caixa | 1000 |
| Unidades por palete | 80000 |
| Profundidade da embalagem externa | 378 mm |
| Largura da embalagem externa | 288 mm |
| Altura da embalagem externa | 180 mm |
| Volume da embalagem externa | 0.0196 cbm |
| Peso do produto | 0.002 kg |
| Peso da embalagem externa | 2.35 kg |
| EAN da caixa interna | 4038917268264 |
| EAN da embalagem externa | 4038917282697 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT