



Product Specification
**Microvette® CB 300 EDTA K2, 300 µl, tampa
violeta, tampa de pressão, fundo plano**

Page 1



Product description

Order number	16.444.100
Descrição do artigo	Microvette® CB 300 EDTA K2, Medicina veterinária, preparo: K2 EDTA, volume nominal: 300 µl, (CxØ) com tampa: 46,5 x 10,8 mm, tampa violeta, tampa de pressão, fundo plano, recipiente interno: em forma de trompeta, 100 unid./pacote

Product characteristics

Tipo de retirada	capilar
Tipo de preparo	K2 EDTA
Código de cor	ISO
Tipo de fechamento	tampa de pressão
Campo de aplicação	Medicina veterinária
Forma do fundo	fundo plano
Forma do recipiente interior	em forma de trompeta

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**Microvette® CB 300 EDTA K2, 300 µl, tampa
violeta, tampa de pressão, fundo plano**

Page 2

Size

Volume nominal	300 µl
Diâmetro	10.8 mm
Comprimento com tampa	46.5 mm
Comprimento sem tampa	43.7 mm

Material & colours

Material do produto	Polipropileno (PP)
Cor do produto	transparente
Material de fechamento	Low Density Polyethylene (LD-PE)
Cor da tampa	violeta

Purity & certification

Categoria do produto	Não é um dispositivo médico DIV
Em lotes	sim

Packaging

Quantidade mínima do pedido	1000
Tipo de embalagem interior	pacote
Unidades por caixa interna	100
Unidades por caixa	1000
Unidades por palete	80000
Profundidade da embalagem externa	378 mm
Largura da embalagem externa	288 mm
Altura da embalagem externa	180 mm
Volume da embalagem externa	0.0196 cbm
Peso do produto	0.002 kg
Peso da embalagem externa	2.39 kg
EAN da caixa interna	4038917268356
EAN da embalagem externa	4038917282772

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT