



Product Specification

Tubo de rosca, 2,5 ml, (CxØ): 75 x 13 mm, fundo falso cônico, fundo do tubo arredondado, PP, sem tampa, 1.000 unid./pacote

Page 1



Product description

| | |
|---------------------|--|
| Order number | 60.614.010 |
| Descrição do artigo | Tubo de rosca, volume de trabalho: 2,5 ml, (CxØ): 75 x 13 mm, fundo falso cônico, fundo do tubo arredondado, transparente, material: PP, com escala, sem tampa, 1.000 unid./pacote, 2.000 unidades/caixa |

Product characteristics

| | |
|--------------------|---|
| Graduação | sim |
| Forma do fundo | fundo falso cônico, fundo do tubo arredondado |
| Volume de trabalho | 2,5 ml |

Size

| | |
|------------------------|-------|
| Diâmetro | 13 mm |
| Altura do produto | 75 mm |
| Comprimento do produto | 75 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Tubo de rosca, 2,5 ml, (CxØ): 75 x 13 mm, fundo falso cônico, fundo do tubo arredondado, PP, sem tampa, 1.000 unid./pacote

Page 2

Material & colours

| | |
|---------------------|--------------------|
| Material do produto | Polipropileno (PP) |
| Cor do produto | transparente |
| Material do tubo | Polipropileno (PP) |

Purity & certification

| | |
|----------------------|--------------------------|
| Categoria do produto | Diagnóstico in vitro, CE |
| Certificação CE | CE |
| Em lotes | sim |

Packaging

| | |
|-----------------------------------|---------------|
| Quantidade mínima do pedido | 2000 |
| Tipo de embalagem interior | pacote |
| Unidades por caixa interna | 1000 |
| Unidades por caixa | 2000 |
| Unidades por palete | 48000 |
| Profundidade da embalagem externa | 458 mm |
| Largura da embalagem externa | 408 mm |
| Altura da embalagem externa | 260 mm |
| Volume da embalagem externa | 0.0486 cbm |
| Peso do produto | 0.0024 kg |
| Peso da embalagem externa | 5.45 kg |
| EAN da caixa interna | 4038917353182 |
| EAN da embalagem externa | 4038917353175 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT