



Product Specification

Tubo de rosca, 5 ml, (CxØ): 92 x 15,3 mm, fundo falso cônico, fundo do tubo arredondado, PP, tampa montada, 100 unid./pacote

Page 1



Product description

Order number	62.612
Descrição do artigo	Tubo de rosca, volume de trabalho: 5 ml, (CxØ): 92 x 15,3 mm, fundo falso cônico, fundo do tubo arredondado, transparente, material: PP, com impressão, rótulo/impressão: branca, com escala, tampa montada, 100 unid./pacote, 1.000 unidades/caixa

Product characteristics

Rótulo/impressão	com impressão
Cor impressão/rótulo	branca
Graduação	sim
Forma do fundo	fundo falso cônico, fundo do tubo arredondado
Volume de trabalho	5 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Tubo de rosca, 5 ml, (CxØ): 92 x 15,3 mm, fundo falso cônico, fundo do tubo arredondado, PP, tampa montada, 100 unid./pacote

Page 2

Size

Diâmetro	15.3 mm
Comprimento do produto	92 mm
Comprimento sem tampa	92 mm

Material & colours

Material do produto	Polipropileno (PP)
Cor do produto	transparente

Purity & certification

Corresponde às exigências	IATA, ADR
Categoria do produto	Diagnóstico in vitro, CE
Certificação CE	CE
Em lotes	sim

Packaging

Quantidade mínima do pedido	1000
Tipo de embalagem interior	pacote
Unidades por caixa interna	100
Unidades por caixa	1000
Unidades por palete	36000
Profundidade da embalagem externa	398 mm
Largura da embalagem externa	398 mm
Altura da embalagem externa	300 mm
Volume da embalagem externa	0.0475 cbm
Peso do produto	0.005 kg
Peso da embalagem externa	5.69 kg
EAN da caixa interna	4038917082884
EAN da embalagem externa	4038917001908

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT