



Product Specification

Kit BKF Safety-Multifly® 21G, 21G x 3/4", verde,
comprimento do tubo flexível: 200 mm, 1
unid./blister

Page 1



Product description

Order number

51.1638.940

Descrição do artigo

Kit BKF Safety-Multifly® 21G, cânula (ØxC): 21G x 3/4", cor de codificação: verde, comprimento do tubo flexível: 200 mm, com kit de adaptadores de hemocultura, componentes do kit: 1 agulha de Safety-Multifly® 21G com tubo flexível longo, isento de DEHP (REF 85.1638.235) e 1 adaptador universal para hemocultura (REF 14.1207), 1 unid./blister, estéril, 1 unid./caixa de papelão

Product characteristics

Versão

kit de adaptadores de hemocultura

Cor de codificação

verde

Produto Safety

sim

Conexão

frasco de cultura de sangue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Kit BKF Safety-Multifly® 21G, 21G x 3/4", verde,
comprimento do tubo flexível: 200 mm, 1
unid./blister**

Page 2

Size

Diâmetro da cânula (externo)	21 G / 0,8 mm
Diâmetro do tubo flexível	2.3 mm
Comprimento da cânula	19 mm
Comprimento do tubo flexível	200 mm

Purity & certification

Categoria do produto	Produto médico
Padrão de pureza	estéril
Em lotes	sim

Packaging

Quantidade mínima do pedido	50000
Tipo de embalagem interior	blister
Unidades por caixa interna	1
Unidades por caixa	50
Unidades por palete	3200
Profundidade da embalagem externa	398 mm
Largura da embalagem externa	298 mm
Altura da embalagem externa	215 mm
Volume da embalagem externa	0.0255 cbm
Peso do produto	0.0082 kg
Peso da embalagem externa	1.09 kg
EAN da caixa interna	4038917408738
EAN da embalagem externa	4038917408721

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT