



## Product Specification

**S-Monovette® RNA Exact, máx. 2,4 ml, tampa verde, (CxØ): 100 x 15 mm, com etiqueta de papel**

**Page 1**



### Product description

Order number	01.2048.001
Descrição do artigo	S-Monovette® RNA Exact, preparo: Estabilizador de RNA, máx. 2,4 ml, tampa de rosca com membrana, tampa verde, (CxØ) sem tampa: 100 x 15 mm, com etiqueta de papel, rótulo/impressão: verde, 20 unid./caixa de papelão, estéril

### Product characteristics

Tipo de retirada	venoso
Tipo de preparo	Estabilizador de RNA
Rótulo/impressão	com etiqueta de papel
Cor impressão/rótulo	verde
Tampa	tampa de rosca com membrana
Tipo de fechamento	tampa de rosca
Hermético	sim
Conexão	tampa de rosca com membrana

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® RNA Exact, máx. 2,4 ml, tampa verde, (CxØ): 100 x 15 mm, com etiqueta de papel**

Page 2

### Size

Volume máximo	2.4 ml
Diâmetro externo do fecho	18 mm
Diâmetro	15 mm
Comprimento com tampa	116 mm
Comprimento sem tampa	100 mm
Comprimento com tampa e barra	136 mm

### Material & colours

Material do produto	Polipropileno (PP)
Cor do produto	transparente
Material de fechamento	High Density Polyethylene (HD-PE)
Cor da tampa	verde
Material da membrana	borracha
Material do tubo	Polipropileno (PP)
Material do êmbolo	Poliétileno de alta densidade (HD-PE)
Material da haste do êmbolo	Poliestireno (PS)

### Purity & certification

Corresponde às exigências	IATA, ADR
Categoria do produto	Diagnóstico in vitro, CE
Certificação CE	CE
Padrão de pureza	estéril
Esterilização	Irradiação de elétrons
Em lotes	sim

### Packaging

Quantidade mínima do pedido	20
Tipo de embalagem interior	caixa de papelão
Unidades por caixa interna	20
Unidades por caixa	80
Unidades por palete	8640
Profundidade da caixa interna	184 mm
Largura da caixa interna	141 mm
Altura da caixa interna	83 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® RNA Exact, máx. 2,4 ml, tampa verde, (CxØ): 100 x 15 mm, com etiqueta de papel**

**Page 3**

Profundidade da embalagem externa	402 mm
Largura da embalagem externa	297 mm
Altura da embalagem externa	100 mm
Volume da embalagem externa	0.0119 cbm
Peso do produto	0.0149 kg
Peso da embalagem externa	1.86 kg
EAN da caixa interna	4038917550475
EAN da embalagem externa	4038917550468

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**