



Product Specification

S-Monovette® Soro com Gel CAT, 9 ml, tampa marrom, (CxØ): 92 x 16 mm, com etiqueta de papel

Page 2

Tempo de centrifugação	10 Min
Hermético	sim
Conexão	tampa de rosca com membrana

Size

Volume nominal	9 ml
Diâmetro externo do fecho	19 mm
Diâmetro	16 mm
Comprimento com tampa	108 mm
Comprimento sem tampa	92 mm
Comprimento com tampa e barra	129 mm

Material & colours

Cor do produto	transparente
Material de fechamento	High Density Polyethylene (HD-PE)
Cor da tampa	marrom
Material da membrana	borracha
Material do tubo	Polipropileno (PP)
Material do êmbolo	Poliétileno de alta densidade (HD-PE)
Material da haste do êmbolo	Poliestireno (PS)

Purity & certification

Corresponde às exigências	IATA, ADR
Categoria do produto	Diagnóstico in vitro, CE
Certificação CE	CE
Padrão de pureza	estéril
Esterilização	Irradiação de elétrons
Em lotes	sim

Packaging

Quantidade mínima do pedido	500
Tipo de embalagem interior	caixa de papelão
Unidades por caixa interna	50
Unidades por caixa	500
Unidades por palete	24000

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Soro com Gel CAT, 9 ml, tampa marrom, (CxØ): 92 x 16 mm, com etiqueta de papel

Page 3

Profundidade da caixa interna	204 mm
Largura da caixa interna	134 mm
Altura da caixa interna	88 mm
Profundidade da embalagem externa	648 mm
Largura da embalagem externa	303 mm
Altura da embalagem externa	150 mm
Volume da embalagem externa	0.0295 cbm
Peso do produto	0.0093 kg
Peso da embalagem externa	5.62 kg
EAN da caixa interna	4038917128087
EAN da embalagem externa	4038917021579

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT