





## Product Specification

**S-Monovette® Soro com Gel CAT, 2,6 ml, tampa marrom, (CxØ): 65 x 13 mm, com etiqueta de papel**

Page 2

Tempo de centrifugação	10 Min
Hermético	sim
Conexão	tampa de rosca com membrana

### Size

Volume nominal	2.6 ml
Diâmetro externo do fecho	16 mm
Diâmetro	13 mm
Comprimento com tampa	81 mm
Comprimento sem tampa	65 mm
Comprimento com tampa e barra	102 mm

### Material & colours

Cor do produto	transparente
Material de fechamento	High Density Polyethylene (HD-PE)
Cor da tampa	marrom
Material da membrana	borracha
Material do tubo	Polipropileno (PP)
Material do êmbolo	Poliétileno de alta densidade (HD-PE)
Material da haste do êmbolo	Poliestireno (PS)

### Purity & certification

Corresponde às exigências	IATA, ADR
Categoria do produto	Diagnóstico in vitro, CE
Certificação CE	CE
Padrão de pureza	estéril
Esterilização	Irradiação de elétrons
Em lotes	sim

### Packaging

Quantidade mínima do pedido	500
Tipo de embalagem interior	caixa de papelão
Unidades por caixa interna	50
Unidades por caixa	500
Unidades por palete	42000

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Soro com Gel CAT, 2,6 ml, tampa marrom, (CxØ): 65 x 13 mm, com etiqueta de papel**

**Page 3**

Profundidade da caixa interna	171 mm
Largura da caixa interna	105 mm
Altura da caixa interna	73 mm
Profundidade da embalagem externa	388 mm
Largura da embalagem externa	358 mm
Altura da embalagem externa	125 mm
Volume da embalagem externa	0.0174 cbm
Peso do produto	0.0049 kg
Peso da embalagem externa	3.15 kg
EAN da caixa interna	4038917134248
EAN da embalagem externa	4038917025553

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**