



Product description

Order number	04.1914.100
Descrição do artigo	S-Monovette® EDTA K3, preparo: K3 EDTA, 3,4 ml, tampa de rosca com membrana, tampa violeta, código de cor ISO, (CxØ) sem tampa: 65 x 13 mm, com etiqueta de papel, rótulo/impressão: violeta, 50 unid./caixa de papelão, estéril

Product characteristics

Tipo de retirada	venoso
Tipo de preparo	K3 EDTA
Concentração do preparo	1,6 mg/ml sangue
Código de cor	ISO
Rótulo/impressão	com etiqueta de papel
Cor impressão/rótulo	violeta
Graduação	marcação nominal
Tampa	tampa de rosca com membrana
Tipo de fechamento	tampa de rosca
Hermético	sim
Conexão	tampa de rosca com membrana

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**S-Monovette® EDTA K3, 3,4 ml, tampa violeta,
(CxØ): 65 x 13 mm, com etiqueta de papel**

Page 2

Size

Volume nominal	3.4 ml
Diâmetro externo do fecho	16 mm
Diâmetro	13 mm
Comprimento com tampa	81 mm
Comprimento sem tampa	65 mm
Comprimento com tampa e barra	102 mm

Material & colours

Cor do produto	transparente
Material de fechamento	High Density Polyethylene (HD-PE)
Cor da tampa	violeta
Material da membrana	borracha
Material do tubo	Polipropileno (PP)
Material do êmbolo	Polietileno de alta densidade (HD-PE)
Material da haste do êmbolo	Poliestireno (PS)

Purity & certification

Corresponde às exigências	IATA, ADR
Categoria do produto	Diagnóstico in vitro, CE
Certificação CE	CE
Padrão de pureza	estéril
Esterilização	Irradiação de eletrões
Em lotes	sim

Packaging

Quantidade mínima do pedido	500
Tipo de embalagem interior	caixa de papelão
Unidades por caixa interna	50
Unidades por caixa	500
Unidades por palete	42000
Profundidade da caixa interna	171 mm
Largura da caixa interna	105 mm
Altura da caixa interna	73 mm
Profundidade da embalagem externa	388 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® EDTA K3, 3,4 ml, tampa violeta,
(CxØ): 65 x 13 mm, com etiqueta de papel

Page 3

Largura da embalagem externa	358 mm
Altura da embalagem externa	125 mm
Volume da embalagem externa	0.0174 cbm
Peso do produto	0.0042 kg
Peso da embalagem externa	2.78 kg
EAN da caixa interna	4038917134361
EAN da embalagem externa	4038917041874

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT