



Product Specification

Tubo de amostra, EDTA K3, 5 ml, tampa vermelha, (CxØ): 75 x 13 mm, com impressão

Page 1



Product description

Order number	32.332
Descrição do artigo	Tubo de amostra, volume nominal: 5 ml, tampa vermelha, preparo: K3 EDTA, tampa de pressão, (CxØ): 75 x 13 mm, fundo redondo, com impressão, 100 unid./pacote

Product characteristics

Tipo de preparo	K3 EDTA
Concentração do preparo	1,6 mg/ml sangue
Código de cor	EU
Rótulo/impressão	com impressão
Cor impressão/rótulo	vermelha
Graduação	sim
Tampa	Tampa com saia
Tipo de fechamento	tampa de pressão
Forma do fundo	fundo redondo

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Tubo de amostra, EDTA K3, 5 ml, tampa vermelha, (CxØ): 75 x 13 mm, com impressão

Page 2

Size

Volume nominal	5 ml
Diâmetro	13 mm
Comprimento do produto	75 mm
Comprimento sem tampa	75 mm

Material & colours

Cor do produto	transparente
Material de fechamento	Low Density Polyethylene (LD-PE)
Cor da tampa	vermelha
Material do tubo	Polipropileno (PP)

Purity & certification

Categoria do produto	Diagnóstico in vitro, CE
Certificação CE	CE
Em lotes	sim

Packaging

Quantidade mínima do pedido	1000
Tipo de embalagem interior	pacote
Unidades por caixa interna	100
Unidades por caixa	1000
Unidades por palete	48000
Profundidade da embalagem externa	388 mm
Largura da embalagem externa	378 mm
Altura da embalagem externa	215 mm
Volume da embalagem externa	0.0315 cbm
Peso do produto	0.0036 kg
Peso da embalagem externa	4.12 kg
EAN da caixa interna	4038917081184
EAN da embalagem externa	4038917000321

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT