



## Product Specification

**Tubo de amostra, Heparina de lítio LH, 10 ml,  
tampa laranja, (CxØ): 101 x 16,5 mm, com etiqueta  
de papel**

**Page 1**



### Product description

|                     |  |
|---------------------|--|
| Order number        | 46.363.001   |
| Descrição do artigo | Tubo de amostra, volume nominal: 10 ml, tampa laranja, preparo: Heparina, tampa de rosca, (CxØ): 101 x 16,5 mm, fundo plano, com etiqueta de papel, 100 unid./pacote |

### Product characteristics

|                           |                       |
|---------------------------|-----------------------|
| Tipo de preparo           | Heparina              |
| Concentração do preparo   | 16 UI/ml sangue       |
| Código de cor             | EU                    |
| Rótulo/impressão          | com etiqueta de papel |
| Cor impressão/rótulo      | laranja               |
| Graduação                 | marcação nominal      |
| Tipo de fechamento        | tampa de rosca        |
| Forma do fundo            | fundo plano           |
| Centrifugação recomendada | 2000 x g              |
| Tempo de centrifugação    | 10 Min                |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Tube de amostra, Heparina de lítio LH, 10 ml,  
tampa laranja, (CxØ): 101 x 16,5 mm, com etiqueta  
de papel**

**Page 2**

### Size

|                        |         |
|------------------------|---------|
| Volume nominal         | 10 ml   |
| Diâmetro               | 16.5 mm |
| Comprimento do produto | 101 mm  |
| Comprimento sem tampa  | 101 mm  |

### Material & colours

|                        |                                  |
|------------------------|----------------------------------|
| Cor do produto         | transparente                     |
| Material de fechamento | Low Density Polyethylene (LD-PE) |
| Cor da tampa           | laranja                          |
| Material do tubo       | Polipropileno (PP)               |

### Purity & certification

|                           |                          |
|---------------------------|--------------------------|
| Corresponde às exigências | IATA, ADR                |
| Categoria do produto      | Diagnóstico in vitro, CE |
| Certificação CE           | CE                       |
| Em lotes                  | sim                      |

### Packaging

|                                      |               |
|--------------------------------------|---------------|
| Quantidade mínima do pedido          | 1000          |
| Tipo de embalagem interior           | pacote        |
| Unidades por caixa interna           | 100           |
| Unidades por caixa                   | 1000          |
| Unidades por palete                  | 30000         |
| Profundidade da embalagem<br>externa | 398 mm        |
| Largura da embalagem<br>externa      | 398 mm        |
| Altura da embalagem externa          | 340 mm        |
| Volume da embalagem<br>externa       | 0.0539 cbm    |
| Peso do produto                      | 0.0076 kg     |
| Peso da embalagem externa            | 8.32 kg       |
| EAN da caixa interna                 | 4038917105736 |
| EAN da embalagem externa             | 4038917010818 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**