



## Product Specification

Recipiente de proteção, forma: quadrangular,  
comprimento: 117 mm, Ø da abertura: 28 mm,  
tampa incluída

Page 1



### Product description

Order number	78.572
Descrição do artigo	Recipiente de proteção, quadrangular, forma: quadrangular, (CxL): 117 x 35 mm, Ø da abertura: 28 mm, tampa incluída, tampa: natural, 50 unid./pacote

### Product characteristics

Tampa	tampa redonda
Tipo de fechamento	tampa de rosca
Campo de aplicação	Envio de amostras
Exemplo de utilização	Embalagem secundária em conformidade com P 650
Compressa de absorção	não
Forma da abertura	redondo
Hermético	sim

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-09-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Recipiente de proteção, forma: quadrangular,  
comprimento: 117 mm, Ø da abertura: 28 mm,  
tampa incluída**

**Page 2**

### Size

Diâmetro da abertura	28 mm
Largura do produto	35 mm
Comprimento do produto	117 mm

### Material & colours

Material do produto	Polietileno de alta densidade (HD-PE)
Cor da tampa	natural

### Purity & certification

Corresponde às exigências	P 650, ADR
Categoria do produto	Não é um dispositivo médico   DIV
Em lotes	sim

### Packaging

Quantidade mínima do pedido	250
Tipo de embalagem interior	pacote
Unidades por caixa interna	50
Unidades por caixa	250
Unidades por palete	6000
Profundidade da embalagem externa	398 mm
Largura da embalagem externa	398 mm
Altura da embalagem externa	415 mm
Volume da embalagem externa	0.0657 cbm
Peso do produto	0.0113 kg
Peso da embalagem externa	3.57 kg
EAN da caixa interna	4038917082587
EAN da embalagem externa	4038917001663

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-09-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**