



### **Product description**

Order number	75.1437.100
Descrição do artigo	Copos para urina, volume de trabalho máx.: 75 ml, (ØxA): 65 x 48 mm, transparente, graduado, material: PP

### **Product characteristics**

Graduação	sim
Volume de trabalho	75 ml

### **Size**

Volume nominal	75 ml
Volume máximo	75 ml
Diâmetro	65 mm
Altura do produto	48 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Copos para urina, 75 ml, (ØxA): 65 x 48 mm, PP,  
transparente**

**Page 2**

**Material & colours**

Material do produto Polipropileno (PP)  
Cor do produto transparente

**Purity & certification**

Categoria do produto Não é um dispositivo médico | DIV  
Em lotes sim

**Packaging**

Quantidade mínima do pedido 1000  
Unidades por caixa interna 25  
Unidades por caixa 1000  
Unidades por palete 28000  
Volume da embalagem 0.065 cbm  
externa  
Peso do produto 0.0092 kg  
Peso da embalagem externa 9.65 kg  
EAN da caixa interna 4038917129138  
EAN da embalagem externa 4038917074711

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-06-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**