



Product Specification

Kit para urina de jato médio, volume máx.: 70 ml,
(CxØ): 115 x 95 mm, para recipiente de 70 ml,
graduado, PP, translúcido, com etiqueta de papel

Page 1



Product description

Order number	51.9923.812
Descrição do artigo	Kit para urina de jato médio, volume máx.: 70 ml, (CxØ): 115 x 95 mm, Ø da abertura: 48,45 mm, para recipiente de 70 ml, translúcido, graduado, material: PP, tampa de rosca, tampa incluída, com etiqueta de papel, estéril, 1 unid./blister

Product characteristics

Rótulo/impressão	com etiqueta de papel
Graduação	sim
Tampa	tampa de rosca padrão
Tipo de fechamento	tampa de rosca
Campo de aplicação	coleta e armazenamento de urina
Exemplo de utilização	Colheita de urina através de um funil para dentro de um recipiente de 70 ml
Indicado para	recipiente de 70 ml
Forma do fundo	fundo plano
Volume de trabalho	70 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Kit para urina de jato médio, volume máx.: 70 ml,
(CxØ): 115 x 95 mm, para recipiente de 70 ml,
graduado, PP, translúcido, com etiqueta de papel**

Page 2

Size

Volume nominal	70 ml
Volume máximo	70 ml
Diâmetro	95 mm
Diâmetro da abertura	48,45 mm
Altura do produto	115 mm
Comprimento do produto	115 mm
Comprimento sem tampa	110 mm

Material & colours

Material do produto	Polipropileno (PP)
Cor do produto	translúcido
Material de fechamento	High Density Polyethylene (HD-PE)

Purity & certification

Isento de látex	sim
Padrão de pureza	estéril
Isenção de pirogênio	não
Ficha de dados de segurança	não

Packaging

Quantidade mínima do pedido	36
Tipo de embalagem interior	blister
Unidades por caixa interna	1
Unidades por caixa	36
Unidades por palete	1440
Profundidade da embalagem externa	604 mm
Largura da embalagem externa	264 mm
Altura da embalagem externa	373 mm
Volume da embalagem externa	0.06 cbm
Peso da embalagem externa	4.4 kg
EAN da caixa interna	4038917154529
EAN da embalagem externa	4038917036733

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT