



### **Product description**

Order number	92.180.614
Descrição do artigo	Rotor em bloco, para 7 tubos até 28 mm Ø, para SARMIX® M 2000, 1 unid./caixa de papelão

### **Product characteristics**

Indicado para	SARMIX® M 2000
Capacidade de armazenamento	7

### **Purity & certification**

Categoria do produto	Não é um dispositivo médico   DIV
Em lotes	não

### **Packaging**

Quantidade mínima do pedido	1
Tipo de embalagem interior	caixa de papelão
Unidades por caixa interna	1
Unidades por caixa	10
Unidades por palete	10

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Rotor em bloco, para 7 tubos até 28 mm Ø, para**  
**SARMIX® M 2000**

**Page 2**

Volume da embalagem externa	0.0122 cbm
Peso do produto	0.2009 kg
Peso da embalagem externa	2.21 kg
EAN da caixa interna	4038917102025
EAN da embalagem externa	4038917075787

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**