



Product description

Order number	92.180.616
Descrição do artigo	Rotor em bloco, para 14 tubos até 28 mm Ø, para SARMIX® M 2000, 1 unid./caixa de papelão

Product characteristics

Indicado para	SARMIX® M 2000
Capacidade de armazenamento	14

Purity & certification

Categoria do produto	Não é um dispositivo médico DIV
Em lotes	não

Packaging

Quantidade mínima do pedido	1
Tipo de embalagem interior	caixa de papelão
Unidades por caixa interna	1
Unidades por caixa	10
Unidades por palete	10

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Rotor em bloco, para 14 tubos até 28 mm Ø, para
SARMIX® M 2000

Page 2

Volume da embalagem externa	0.0122 cbm
Peso do produto	0.18 kg
Peso da embalagem externa	3.01 kg
EAN da caixa interna	4038917160858
EAN da embalagem externa	4038917160841

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT