




Product description

Order number 92.181.610
Descrição do artigo Rotor, para 36 recipientes, para tubos Ø 10,8 mm

Product characteristics

Indicado para tubos Ø 10,8 mm
Capacidade de armazenamento 36

Purity & certification

Categoria do produto Não é um dispositivo médico | DIV
Em lotes sim

Packaging

Quantidade mínima do pedido 1
Unidades por caixa interna 1
Unidades por caixa 1
Unidades por palete 40
Volume da embalagem 0.0056 cbm
externa
Peso do produto 1.3595 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Rotor, para 36 recipientes, para tubos Ø 10,8 mm

Page 2

Peso da embalagem externa 1.36 kg
EAN da caixa interna 4038917156233
EAN da embalagem externa 4038917156226

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT