



### Product description

Order number 92.181.612  
Descrição do artigo Rotor, para 6 recipientes, para tubos Ø 28 mm

### Product characteristics

Indicado para tubos Ø 28 mm  
Capacidade de armazenamento 6

### Purity & certification

Categoria do produto Não é um dispositivo médico | DIV  
Em lotes sim

### Packaging

Quantidade mínima do pedido 1  
Unidades por caixa interna 1  
Unidades por caixa 1  
Unidades por palete 40  
Volume da embalagem 0.0053 cbm  
externa  
Peso do produto 1.3595 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Rotor, para 6 recipientes, para tubos Ø 28 mm**

**Page 2**

Peso da embalagem externa 1.36 kg  
EAN da caixa interna 4038917156271  
EAN da embalagem externa 4038917156264

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**