



Product Specification

Cotonete de esfregação para ciência forense,
redondo, no tubo com membrana de ventilação,
ISO 18385, 95 mm, algodão

Page 1



Product description

Order number	80.629.001
Descrição do artigo	Cotonete de esfregação para ciência forense, redondo, no tubo com membrana de ventilação, com etiqueta de papel, tratada com ETO, ISO 18385, Secagem passiva, comprimento da barra: 95 mm, material da barra: madeira, cabeça do toalhete redondo, material do toalhete: algodão, tubo: PP, blister livre de alumínio, 1 unid./blister

Product characteristics

Rótulo/impressão	com etiqueta de papel
Tampa	tampa de rosca padrão
Tipo de fechamento	tampa de rosca
Versão	no tubo com membrana de ventilação
Campo de aplicação	Coleta de amostras forenses

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Cotonete de esfregaço para ciência forense,
redondo, no tubo com membrana de ventilação,
ISO 18385, 95 mm, algodão**

Page 2

Size

Diâmetro	15.35 mm
Largura do produto	18.05 mm
Altura do produto	106.8 mm
Comprimento do produto	107 mm
Comprimento com tampa	97.6 µm
Comprimento sem tampa	106 mm
Comprimento com tampa e barra	106.8 µm
Comprimento da barra	95 mm

Material & colours

Cor do produto	de várias cores
Material de fechamento	High Density Polyethylene (HD-PE)
Cor da tampa	branca
Material do tampão	algodão
Material da barra	madeira
Material do tubo	Polipropileno (PP)

Purity & certification

Corresponde às exigências	ISO 18385
Categoria do produto	Não é um dispositivo médico DIV
Padrão de pureza	tratada com ETO
Esterilização	Esterilização a gás
Em lotes	sim

Packaging

Quantidade mínima do pedido	800
Tipo de embalagem interior	blister
Unidades por caixa interna	100
Unidades por caixa	800
Unidades por palete	20000
Profundidade da caixa interna	211 mm
Largura da caixa interna	151 mm
Altura da caixa interna	171 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Cotonete de esfregação para ciência forense,
redondo, no tubo com membrana de ventilação,
ISO 18385, 95 mm, algodão**

Page 3

Profundidade da embalagem externa	437 mm
Largura da embalagem externa	363 mm
Altura da embalagem externa	335 mm
Volume da embalagem externa	0.0531 cbm
Peso do produto	0.0065 kg
Peso da embalagem externa	7.29 kg
EAN da caixa interna	4038917276887
EAN da embalagem externa	4038917428941

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT