



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 94.6077.136  |
| Descrição do artigo | miniPERM® , Seringa descartável 2 ml, para Biorreator miniPERM® , Luer, não autoclavável, estéril, 1 unid./blister |

### **Product characteristics**

|               |                      |
|---------------|----------------------|
| Indicado para | Biorreator miniPERM® |
| Conexão       | Luer                 |

### **Size**

|               |      |
|---------------|------|
| Volume máximo | 2 ml |
|---------------|------|

### **Material & colours**

|                     |                    |
|---------------------|--------------------|
| Material do produto | Polipropileno (PP) |
|---------------------|--------------------|

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**miniPERM® , Seringa descartável 2 ml, para**  
**Biorreator miniPERM® , Luer**

**Page 2**

**Purity & certification**

|                      |                     |
|----------------------|---------------------|
| Categoria do produto | Produto médico      |
| Certificação CE      | Fornecedor CE       |
| Padrão de pureza     | estéril             |
| Esterilização        | Esterilização a gás |
| Isenção de pirogênio | sim                 |
| Em lotes             | sim                 |

**Packaging**

|                                   |               |
|-----------------------------------|---------------|
| Quantidade mínima do pedido       | 100           |
| Tipo de embalagem interior        | blister       |
| Unidades por caixa interna        | 100           |
| Unidades por caixa                | 100           |
| Unidades por palete               | 32400         |
| Profundidade da embalagem externa | 268 mm        |
| Largura da embalagem externa      | 191 mm        |
| Altura da embalagem externa       | 100 mm        |
| Volume da embalagem externa       | 0.0051 cbm    |
| Peso do produto                   | 0.003 kg      |
| Peso da embalagem externa         | 0.42 kg       |
| EAN da caixa interna              | 4038917255929 |
| EAN da embalagem externa          | 4038917255912 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**