



## Product Specification

**Tubo para a cultura de células, (CxØ): 125 x 16 mm, fundo redondo, com tratamento TC**

**Page 1**



### Product description

Order number	83.9923.945
Descrição do artigo	Tubo para a cultura de células, 15 ml, (CxØ): 125 x 16 mm, fundo redondo, PS, com tratamento TC, para Culturas de suspensões e de monocamadas, estéril, 5 unid./pacote

### Product characteristics

Graduação	não
Tampa	tampa de rosca padrão
Tipo de célula	Culturas de suspensões e de monocamadas
Forma do fundo	fundo redondo
Volume de trabalho	15 ml
Tipo de recipiente	Tubo
Tratamento da superfície	com tratamento TC
Cor de codificação	vermelha

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Tubo para a cultura de células, (CxØ): 125 x 16 mm, fundo redondo, com tratamento TC**

**Page 2**

### Size

Diâmetro	16 mm
Comprimento do produto	125 mm

### Material & colours

Material do produto	Poliestireno (PS)
Cor do produto	transparente
Material de fechamento	Polyethylene (PE)
Cor da tampa	vermelha

### Purity & certification

Categoria do produto	Não é um dispositivo médico   DIV
Padrão de pureza	estéril
Esterilização	Radiação gama
Em lotes	sim

### Packaging

Quantidade mínima do pedido	1000
Tipo de embalagem interior	pacote
Unidades por caixa interna	5
Unidades por caixa	1000
Unidades por palete	21000
Profundidade da embalagem externa	428 mm
Largura da embalagem externa	438 mm
Altura da embalagem externa	372 mm
Volume da embalagem externa	0.0697 cbm
Peso do produto	0.0086 kg
Peso da embalagem externa	9.4 kg
EAN da caixa interna	4038917154635
EAN da embalagem externa	4038917027892

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**