



## Product description

Order number	92.170.390
Descrição do artigo	Maleta de transporte GS 100, para Coletor de germes do ar, 1 unid./caixa de papelão

## Product characteristics

Indicado para	Coletor de germes do ar
---------------	-------------------------

## Purity & certification

Categoria do produto	Não é um dispositivo médico   DIV
Em lotes	não

## Packaging

Quantidade mínima do pedido	1
Tipo de embalagem interior	caixa de papelão
Unidades por caixa interna	1
Unidades por caixa	1
Unidades por palete	1
Volume da embalagem externa	0.0286 cbm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Maleta de transporte GS 100, para Coletor de germes do ar**

**Page 2**

Peso do produto	2.7395 kg
Peso da embalagem externa	2.74 kg
EAN da caixa interna	4038917101295
EAN da embalagem externa	4038917051590

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**