



### **Product description**

Order number	90.189.220
Descrição do artigo	Filtro de reposição, para Automatic-Sarpette® , tamanho do poro: 0,45 µm, 5 unid./pacote

### **Product characteristics**

Indicado para	Automatic-Sarpette®
Tamanho do poro	0,45

### **Purity & certification**

Categoria do produto	Não é um dispositivo médico   DIV
Em lotes	não

### **Packaging**

Quantidade mínima do pedido	5
Tipo de embalagem interior	pacote
Unidades por caixa interna	5
Unidades por caixa	500
Unidades por palete	54000

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Filtro de reposição, para Automatic-Sarpette® ,**  
**tamanho do poro: 0,45 µm**

**Page 2**

Profundidade da embalagem externa	308 mm
Largura da embalagem externa	223 mm
Altura da embalagem externa	195 mm
Volume da embalagem externa	0.0134 cbm
Peso do produto	0.002 kg
Peso da embalagem externa	1.21 kg
EAN da caixa interna	4038917419024
EAN da embalagem externa	4038917419017

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**