



Product description

| | |
|---------------------|---|
| Order number | 86.1176.300 |
| Descrição do artigo | Pipeta de transferência, volume nominal: 2 ml, (CxL): 154 x 11 mm, material: PEBD, transparente, tamanho da gota: 35 µl - 45 µl, 500 unid./caixa de papelão |

Product characteristics

| | |
|-----------|-----|
| Graduação | não |
|-----------|-----|

Size

| | |
|------------------------|---------------|
| Volume nominal | 2 ml |
| Largura do produto | 11 mm |
| Comprimento do produto | 154 mm |
| Tamanho da gota | 35 µl - 45 µl |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Pipeta de transferência, 2 ml, (CxL): 154 x 11 mm,
PEBD, transparente

Page 2

Material & colours

Material do produto Polietileno de baixa densidade (PEBD)
Cor do produto transparente

Purity & certification

Categoria do produto Não é um dispositivo médico | DIV
Em lotes sim

Packaging

Quantidade mínima do pedido 5000
Tipo de embalagem interior caixa de papelão
Unidades por caixa interna 500
Unidades por caixa 5000
Unidades por palete 84000
Volume da embalagem 0.0994 cbm
externa
Peso do produto 0.0015 kg
Peso da embalagem externa 8.64 kg
EAN da caixa interna 4038917378567
EAN da embalagem externa 4038917378550

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT