



Product Specification
DishRack, altura: 370 mm, amarela, para 88
placas de Petri Ø 92 mm

Page 1



Product description

Order number 93.1646.002
Descrição do artigo DishRack, altura: 370 mm, material: PC, amarela, para 88
placas de Petri Ø 92 mm, 1 unid./caixa de papelão

Product characteristics

Capacidade de armazenamento 88

Size

Altura do produto 370 mm

Material & colours

Material do produto Policarbonato (PC)
Cor do produto amarela

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
DishRack, altura: 370 mm, amarela, para 88
placas de Petri Ø 92 mm

Page 2

Purity & certification

Categoria do produto	Não é um dispositivo médico DIV
Em lotes	não

Packaging

Quantidade mínima do pedido	1
Tipo de embalaje interior	caixa de papelão
Unidades por caixa interna	1
Unidades por caixa	1
Unidades por palete	108
Profundidade da embalagem externa	383 mm
Largura da embalagem externa	193 mm
Altura da embalagem externa	190 mm
Volume da embalagem externa	0.014 cbm
Peso do produto	1.086 kg
Peso da embalagem externa	1.35 kg
EAN da caixa interna	4038917132565
EAN da embalagem externa	4038917030472

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT