



### **Product description**

Order number	65.716.008
Descrição do artigo	Tampa de rosca, violeta, adequado para microtubo com tampa de rosca, 500 unid./pacote

### **Product characteristics**

Tampa	tampa plana
Tipo de fechamento	tampa de rosca
Indicado para	microtubo com tampa de rosca

### **Material & colours**

Material do produto	Polipropileno (PP)
Cor do produto	violeta
Material do O-ring	Borracha de etileno-propileno-dieno (EPDM)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tampa de rosca, violeta, adequado para  
microtubo com tampa de rosca**

**Page 2**

**Purity & certification**

Categoria do produto	Diagnóstico in vitro, CE
Certificação CE	CE
Em lotes	sim

**Packaging**

Quantidade mínima do pedido	10000
Tipo de embalagem interior	pacote
Unidades por caixa interna	500
Unidades por caixa	10000
Unidades por palete	660000
Profundidade da embalagem externa	438 mm
Largura da embalagem externa	268 mm
Altura da embalagem externa	155 mm
Volume da embalagem externa	0.0182 cbm
Peso do produto	0.0003 kg
Peso da embalagem externa	3.52 kg
EAN da caixa interna	4038917118613
EAN da embalagem externa	4038917016261

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**