



Product Specification
Recipiente individual PCR, 0,5 ml, PCR
Performance Tested, cores variadas, PP, tampa
plana

Page 1



Product description

Order number 72.735.992
Descrição do artigo Recipiente individual PCR, 0,5 ml, PCR Performance Tested, vermelho/verde/amarelo/azul/violeta/laranja, material: PP, com proteção anticontaminação, tampa plana, por cor 500 unid./pacote com minialça

Product characteristics

Graduação sim
Tampa tampa plana
Tipo de fechamento tampa de pressão
Indicada para qPCR não
Forma do fundo cônico
High profile ou Low profile Perfil alto
Quantidade de poços 1
Volume de trabalho 25 µl - 0,5 ml
Centrifugação máxima (RCF) 8000 x g
Temperatura de centrifugação 22
Tempo de centrifugação 5 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Recipiente individual PCR, 0,5 ml, PCR
Performance Tested, cores variadas, PP, tampa
plana

Autoclavável (máx. °C) 121 °C

Size

Volume nominal 0.5 ml
Volume máximo 500 µl
Diâmetro 7.9 mm
Largura do produto 10 mm
Altura do produto 29.5 mm
Comprimento do produto 30.9 mm

Material & colours

Material do produto Polipropileno (PP)
Cor do produto cores variadas

Purity & certification

Categoria do produto Diagnóstico in vitro, CE
Certificação CE CE
Padrão de pureza PCR Performance Tested
Em lotes sim

Packaging

Quantidade mínima do pedido 3000
Tipo de embalagem interior pacote com minialça
Unidades por caixa interna 3000
Unidades por caixa 3000
Unidades por palete 198000
Profundidade da embalagem 428 mm
externa
Largura da embalagem 323 mm
externa
Altura da embalagem externa 155 mm
Volume da embalagem 0.0214 cbm
externa
Peso do produto 0.0004 kg
Peso da embalagem externa 1.67 kg
EAN da caixa interna 4038917120067

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Recipiente individual PCR, 0,5 ml, PCR
Performance Tested, cores variadas, PP, tampa
plana

Page 3

EAN da embalagem externa 4038917017053

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT