




## Product Specification

Caixa de armazenamento, tampa articulada, PP,  
dimensão da grade: 10 x 10, para 100 recipientes

Page 1



### Product description

Order number	93.877
Descrição do artigo	Caixa de armazenamento, tampa articulada, material: PP, transparente, dimensão da grade: 10 x 10, para 100 recipientes, adequado para recipientes com dimensões máx. de 45 x 12 mm, 5 unid./pacote

### Product characteristics

Tipo de fechamento	tampa articulada
Indicado para	recipientes com dimensões máx. de 45 x 12 mm
Capacidade de armazenamento	100
Dimensão da grade	10 x 10
Autoclavável (máx. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-06-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

Caixa de armazenamento, tampa articulada, PP,  
dimensão da grade: 10 x 10, para 100 recipientes

Page 2

### Size

Largura do produto	140 mm
Altura do produto	55 mm
Comprimento do produto	150 mm

### Material & colours

Material do produto	Polipropileno (PP)
Cor do produto	transparente
Material de fechamento	Polypropylene (PP)
Cor da tampa	transparente

### Purity & certification

Categoria do produto	Não é um dispositivo médico   DIV
Em lotes	sim

### Packaging

Quantidade mínima do pedido	20
Tipo de embalagem interior	pacote
Unidades por caixa interna	5
Unidades por caixa	20
Unidades por palete	720
Profundidade da embalagem externa	305 mm
Largura da embalagem externa	279 mm
Altura da embalagem externa	305 mm
Volume da embalagem externa	0.026 cbm
Peso do produto	0.1835 kg
Peso da embalagem externa	4.1 kg
EAN da caixa interna	4038917161152
EAN da embalagem externa	4038917161145

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-06-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**