



Product Specification
Rack, poliestireno expandido, dimensão da grade: 5 x 4

Page 1



Product description

Order number	95.64.259
Descrição do artigo	Rack, material: poliestireno expandido, branca, dimensão da grade: 5 x 4, para 20 recipientes, 8 unid./caixa de papelão

Product characteristics

Capacidade de armazenamento	20
Dimensão da grade	5 x 4

Material & colours

Material do produto	poliestireno expandido
Cor do produto	branca

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Rack, poliestireno expandido, dimensão da grade: 5 x 4

Page 2

Purity & certification

Categoria do produto Não é um dispositivo médico | DIV
Em lotes sim

Packaging

Quantidade mínima do pedido 56
Tipo de embalagem interior caixa de papelão
Unidades por caixa interna 8
Unidades por caixa 56
Unidades por palete 1792
Volume da embalagem externa 0.0559 cbm
Peso do produto 0.0238 kg
Peso da embalagem externa 1.36 kg
EAN da caixa interna 4038917386883
EAN da embalagem externa 4038917386876

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT