



Product Specification

Caixa de armazenamento, tampa de encaixe,
papelão, dimensão da grade: 4 x 4, para 16
recipientes

Page 1



Product description

Order number	95.64.916
Descrição do artigo	Caixa de armazenamento, tampa de encaixe, material: papelão, branca, dimensão da grade: 4 x 4, para 16 recipientes, adequado para Tubos de até no máximo 34 mm de Ø, altura do tubo entre 89-115 mm, 1 unid./caixa de papelão

Product characteristics

Tipo de fechamento	tampa de encaixe
Indicado para	Tubos de até no máximo 34 mm de Ø, altura do tubo entre 89-115 mm
Capacidade de armazenamento	16
Dimensão da grade	4 x 4

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-06-30

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Caixa de armazenamento, tampa de encaixe,
papelão, dimensão da grade: 4 x 4, para 16
recipientes

Page 2

Size

Largura do produto	148 mm
Altura do produto	116 mm
Comprimento do produto	148 mm

Material & colours

Material do produto	papelão
Cor do produto	branca
Material de fechamento	cardboard
Cor da tampa	branca

Purity & certification

Categoria do produto	Não é um dispositivo médico DIV
Em lotes	sim

Packaging

Quantidade mínima do pedido	1
Tipo de embalagem interior	caixa de papelão
Unidades por caixa interna	1
Unidades por caixa	30
Unidades por palete	360
Profundidade da embalagem externa	495 mm
Largura da embalagem externa	419 mm
Altura da embalagem externa	546 mm
Volume da embalagem externa	0.1134 cbm
Peso do produto	0.1337 kg
Peso da embalagem externa	5.15 kg
EAN da caixa interna	4038917095679
EAN da embalagem externa	4038917007832

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-06-30

This document was prepared by EDP
support and is valid without signature.



SARSTEDT