




Product Specification

Caixa de armazenamento, tampa de encaixe,
papalão, dimensão da grade: 7 x 7, para 49
recipientes

Page 1



Product description

Order number	95.64.950
Descrição do artigo	Caixa de armazenamento, tampa de encaixe, material: papalão, branca, dimensão da grade: 7 x 7, para 49 recipientes, adequado para tubos Ø 16-17 mm

Product characteristics

Tipo de fechamento	tampa de encaixe
Indicado para	tubos Ø 16-17 mm
Capacidade de armazenamento	49
Dimensão da grade	7 x 7

Size

Largura do produto	135 mm
Altura do produto	87 mm
Comprimento do produto	135 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP
support and is valid without signature.



SARSTEDT



Product Specification

Caixa de armazenamento, tampa de encaixe,
papelão, dimensão da grade: 7 x 7, para 49
recipientes

Page 2

Material & colours

Material do produto	papelão
Cor do produto	branca
Material de fechamento	cardboard
Cor da tampa	branca

Purity & certification

Categoria do produto	Não é um dispositivo médico DIV
Em lotes	sim

Packaging

Quantidade mínima do pedido	60
Unidades por caixa interna	1
Unidades por caixa	60
Unidades por palete	960
Volume da embalagem externa	0.1134 cbm
Peso do produto	0.1373 kg
Peso da embalagem externa	8.27 kg
EAN da caixa interna	4038917319607
EAN da embalagem externa	4038917319591

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT