



Product description

Order number	58.485.300
Descrição do artigo	Tubo, volume de trabalho: 7 ml, (CxØ): 50 x 16 mm, material: PS, fundo plano, transparente, tampa de pressão, 1.000 unid./pacote

Product characteristics

Graduação	não
Tipo de fechamento	tampa de pressão
Forma do fundo	fundo plano
Volume de trabalho	7 ml

Size

Diâmetro	16 mm
Altura do produto	50 mm
Comprimento do produto	50 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tubo, 7 ml, (CxØ): 50 x 16 mm, PS

Page 2

Material & colours

Material do produto	Poliestireno (PS)
Cor do produto	transparente

Purity & certification

Categoria do produto	Não é um dispositivo médico DIV
Em lotes	sim

Packaging

Quantidade mínima do pedido	1000
Tipo de embalagem interior	pacote
Unidades por caixa interna	1000
Unidades por caixa	2000
Unidades por palete	72000
Volume da embalagem externa	0.0492 cbm
Peso do produto	0.0028 kg
Peso da embalagem externa	6.07 kg
EAN da caixa interna	4038917367462
EAN da embalagem externa	4038917367455

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT