



Product Specification

Tubo de rosca, 2,5 ml, (CxØ): 75 x 13 mm, fundo falso cônico, fundo do tubo arredondado, PP, 442 unid./embalagem empilhável

Page 1



Product description

Order number	60.1575.115
Descrição do artigo	Tubo de rosca, volume de trabalho: 2,5 ml, (CxØ): 75 x 13 mm, fundo falso cônico, fundo do tubo arredondado, transparente, material: PP, 442 unid./embalagem empilhável, 3.536 unidades/caixa

Product characteristics

Graduação	não
Tipo de fechamento	tampa de rosca
Forma do fundo	fundo falso cônico, fundo do tubo arredondado
Volume de trabalho	2,5 ml

Size

Diâmetro	13 mm
Comprimento do produto	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Tubo de rosca, 2,5 ml, (CxØ): 75 x 13 mm, fundo falso cônico, fundo do tubo arredondado, PP, 442 unid./embalagem empilhável

Page 2

Material & colours

Material do produto	Polipropileno (PP)
Cor do produto	transparente
Material do tubo	Polipropileno (PP)

Purity & certification

Categoria do produto	Não é um dispositivo médico DIV
Em lotes	não

Packaging

Quantidade mínima do pedido	3536
Tipo de embalagem interior	embalagem empilhável
Unidades por caixa interna	442
Unidades por caixa	3536
Unidades por palete	106080
Profundidade da caixa interna	381 mm
Largura da caixa interna	77 mm
Altura da caixa interna	176 mm
Profundidade da embalagem externa	398 mm
Largura da embalagem externa	358 mm
Altura da embalagem externa	325 mm
Volume da embalagem externa	0.0463 cbm
Peso do produto	0.0024 kg
Peso da embalagem externa	9.56 kg
EAN da caixa interna	4038917340823
EAN da embalagem externa	4038917340816

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT