



Product Specification

Tubo de rosca, 5 ml, (CxØ): 75 x 13 mm, fundo redondo, PP, sem tampa, 500 unid./embalagem empilhável

Page 1



Product description

Order number	60.504.015
Descrição do artigo	Tubo de rosca, volume de trabalho: 5 ml, (CxØ): 75 x 13 mm, fundo redondo, transparente, material: PP, sem tampa, 500 unid./embalagem empilhável, 2.000 unidades/caixa

Product characteristics

Gradação	não
Tipo de fechamento	tampa de rosca
Forma do fundo	fundo redondo
Volume de trabalho	5 ml
Centrifugação máxima (RCF)	5000 x g

Size

Diâmetro	13 mm
Comprimento do produto	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Tubo de rosca, 5 ml, (CxØ): 75 x 13 mm, fundo redondo, PP, sem tampa, 500 unid./embalagem empilhável

Page 2

Material & colours

Material do produto	Polipropileno (PP)
Cor do produto	transparente
Material do tubo	Polipropileno (PP)

Purity & certification

Corresponde às exigências	IATA
Categoria do produto	Diagnóstico in vitro, CE
Certificação CE	CE
Em lotes	sim

Packaging

Quantidade mínima do pedido	2000
Tipo de embalagem interior	embalagem empilhável
Unidades por caixa interna	500
Unidades por caixa	2000
Unidades por palete	112000
Profundidade da caixa interna	413 mm
Largura da caixa interna	76 mm
Altura da caixa interna	189 mm
Profundidade da embalagem externa	428 mm
Largura da embalagem externa	383 mm
Altura da embalagem externa	175 mm
Volume da embalagem externa	0.0287 cbm
Peso do produto	0.0024 kg
Peso da embalagem externa	5.66 kg
EAN da caixa interna	4038917323383
EAN da embalagem externa	4038917323376

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT