



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 19.930.125   |
| Product description | Blood gas capillaries, preparation: calcium-balanced lithium heparin, nominal volume: 125 µl, (LxØ): 110 x 2.05 mm, 250 piece(s)/container |

### **Product characteristics**

|                           |                                  |
|---------------------------|----------------------------------|
| Type of collection        | capillary                        |
| Type of preparation       | calcium-balanced lithium heparin |
| Preparation concentration | 70 IU/ml blood                   |
| Application               | blood gas analysis               |

### **Size**

|                   |         |
|-------------------|---------|
| Sample volume     | 125 µl  |
| Diameter          | 2.05 mm |
| Length of product | 110 mm  |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material Polyethylene terephthalate (PET)  
Colour of product transparent

**Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Latex-free yes  
Batched yes

**Packaging**

Minimum order qty. 250  
Type of smallest container  
subpackaging  
Piece(s) / inner box 250  
Piece (s) / outer case 16000  
Piece(s) / pallet 1280000  
Depth of case 378 mm  
Width of case 288 mm  
Height of case 180 mm  
Case volume 0.0196 cbm  
Weight of product 0.0003 kg  
Weight of case 6.64 kg  
EAN of inner box 4038917234283  
EAN case 4038917234276

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

