



## Product Specification

Faeces tube, with spoon, screw cap, (LxØ): 101 x 16.5 mm, transparent, sterile

Page 1



### Product description

Order number	80.624
Product description	Faeces tube, with spoon, transparent, screw cap, cap: brown, (LxØ) without cap: 101 x 16.5 mm, sterile

### Product characteristics

Closure type	screw cap
Light protection	no
Design	spoon

### Size

Diameter	16.5 mm
Length of product	102 mm
Length excluding cap	101 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Faeces tube, with spoon, screw cap, (LxØ): 101 x 16.5 mm, transparent, sterile**

**Page 2**

**Material & colours**

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Batched	yes

**Packaging**

Minimum order qty.	250
Piece(s) / inner box	50
Piece (s) / outer case	250
Piece(s) / pallet	10500
Depth of case	388 mm
Width of case	378 mm
Height of case	250 mm
Case volume	0.0367 cbm
Weight of product	0.0169 kg
Weight of case	4.69 kg
EAN of inner box	4038917082990
EAN case	4038917001977

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**