




**Product Specification**  
**Faeces tube, screw cap, (LxØ): 76 x 20 mm,**  
**transparent**

**Page 1**



**Product description**

Order number	80.734.301
Product description	Faeces tube, transparent, screw cap, cap: brown, (LxØ) without cap: 76 x 20 mm, with paper label, 100 piece(s)/bag

**Product characteristics**

Label/ Print	with paper label
Colour of print/label	black
Closure type	screw cap

**Size**

Diameter	20 mm
Length of product	77 mm
Length excluding cap	76 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Faeces tube, screw cap, (LxØ): 76 x 20 mm,**  
**transparent**

**Page 2**

**Material & colours**

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	21000
Depth of case	388 mm
Width of case	378 mm
Height of case	235 mm
Case volume	0.0345 cbm
Weight of product	0.0071 kg
Weight of case	4.14 kg
EAN of inner box	4038917119955
EAN case	4038917016957

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**