



### Product description

Order number	16.444.100
Product description	Microvette® CB 300 EDTA K2E, Veterinary medicine, preparation: K2 EDTA, nominal volume: 300 µl, (LxØ) with cap: 46.5 x 10.8 mm, cap violet, push cap, flat base, inner receptacle: trumpet-shaped, 100 piece(s)/bag

### Product characteristics

Type of collection	capillary
Type of preparation	K2 EDTA
Colour code	ISO
Closure type	push cap
Application	Veterinary medicine
Base shape	flat base
Shape of inner tube	trumpet-shaped

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### Size

Sample volume	300 µl
Diameter	10.8 mm
Length including cap	46.5 mm
Length excluding cap	43.7 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	violet

### Purity & certification

Product category	no medical device   IVD
Batched	yes

### Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	80000
Depth of case	378 mm
Width of case	288 mm
Height of case	180 mm
Case volume	0.0196 cbm
Weight of product	0.002 kg
Weight of case	2.39 kg
EAN of inner box	4038917268356
EAN case	4038917282772

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

