



## Product Specification

**Minivette® POCT Lithium heparin LH, 50 µl,  
plunger green, colour code ISO, 200 piece(s)/bag**

**Page 1**



### Product description

Order number	17.2112.150
Product description	Minivette® POCT Lithium heparin LH, nominal volume: 50 µl, preparation: Heparin, green, colour code ISO, 200 piece(s)/bag, 200 piece(s)/case

### Product characteristics

Type of collection	capillary
Type of preparation	Heparin
Colour code	ISO
Application	capillary blood collection
Coding colour	green

### Size

Sample volume	50 µl
---------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Minivette® POCT Lithium heparin LH, 50 µl,  
plunger green, colour code ISO, 200 piece(s)/bag**

**Page 2**

**Material & colours**

Product material Cellulose acetate (CA)  
Colour of product transparent

**Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Latex-free yes  
Batched yes

**Packaging**

Minimum order qty. 2000  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 200  
Piece (s) / outer case 2000  
Piece(s) / pallet 144000  
Depth of box 140 mm  
Width of box 99 mm  
Height of box 117 mm  
Depth of case 513 mm  
Width of case 293 mm  
Height of case 140 mm  
Case volume 0.021 cbm  
Weight of product 0.001 kg  
Weight of case 2.68 kg  
EAN of inner box 4038917320023  
EAN case 4038917320016

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**