



Product Specification

**Minivette® POCT Lithium heparin LH, 200 µl,
plunger green, colour code ISO, 150 piece(s)/bag**

Page 1



Product description

| | |
|---------------------|---|
| Order number | 17.2112.201 |
| Product description | Minivette® POCT Lithium heparin LH, nominal volume: 200 µl, preparation: Heparin, green, colour code ISO, 150 piece(s)/bag, 150 piece(s)/case |

Product characteristics

| | |
|---------------------|----------------------------|
| Type of collection | capillary |
| Type of preparation | Heparin |
| Colour code | ISO |
| Application | capillary blood collection |
| Coding colour | green |

Size

| | |
|---------------|--------|
| Sample volume | 200 µl |
|---------------|--------|

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Material & colours

| | |
|-------------------|------------------------|
| Product material | Cellulose acetate (CA) |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Latex-free | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 150 |
| Piece (s) / outer case | 1500 |
| Piece(s) / pallet | 108000 |
| Depth of box | 140 mm |
| Width of box | 99 mm |
| Height of box | 117 mm |
| Depth of case | 513 mm |
| Width of case | 293 mm |
| Height of case | 140 mm |
| Case volume | 0.021 cbm |
| Weight of product | 0.0011 kg |
| Weight of case | 2.37 kg |
| EAN of inner box | 4038917339513 |
| EAN case | 4038917339506 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

