



### Product description

Order number	17.2113.120
Product description	Minivette® POCT EDTA K3E, nominal volume: 20 µl, preparation: K3 EDTA, violet, colour code ISO, 200 piece(s)/bag, 200 piece(s)/case

### Product characteristics

Type of collection	capillary
Type of preparation	K3 EDTA
Colour code	ISO
Application	capillary blood collection
Coding colour	violet

### Size

Sample volume	20 µl
---------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Cellulose acetate (CA)
Colour of product	transparent

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Latex-free	yes
Batched	yes

### **Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	200
Piece (s) / outer case	2000
Piece(s) / pallet	144000
Depth of box	140 mm
Width of box	99 mm
Height of box	117 mm
Depth of case	513 mm
Width of case	293 mm
Height of case	140 mm
Case volume	0.021 cbm
Weight of product	0.0009 kg
Weight of case	2.53 kg
EAN of inner box	4038917320054
EAN case	4038917320047

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

