



**Product Specification**  
**End-to-end capillary, Ammonium heparin AH, 10  
µl, Ø: 1.25 mm, glass**

**Page 1**

**Product description**

Order number 19.414.001  
Product description End-to-end capillary, Ammonium heparin AH, capillary blood collection, preparation: Ammonium heparin, nominal volume: 10 µl, Ø: 1.25 mm, glass, 100 piece(s)/container

**Product characteristics**

Type of collection capillary  
Type of preparation Ammonium heparin  
Application capillary blood collection  
Volume of work 10 µl

**Size**

Sample volume 10 µl  
Blood volume 10 µl  
Diameter 1.25 mm  
Length of product 29 mm

**Material & colours**

Product material glass  
Colour of product transparent

**Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE supplier  
Batched yes

**Packaging**

Minimum order qty. 1000  
Type of smallest container  
subpackaging  
Piece(s) / inner box 1000  
Piece (s) / outer case 50000  
Piece(s) / pallet 2400000  
Depth of box 117 mm  
Width of box 50 mm  
Height of box 40 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**End-to-end capillary, Ammonium heparin AH, 10  
µl, Ø: 1.25 mm, glass**

Depth of case	388 mm
Width of case	378 mm
Height of case	215 mm
Case volume	0.0315 cbm
Weight of product	0.0001 kg
Weight of case	8.65 kg
EAN of inner box	5706860179225
EAN case	5706860179232

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

