



## Product Specification

Sample tube, Roche HbA1c haemolysis solution,  
1,000 µl, cap black, (LxØ): 82 x 11.5 mm, with  
paper label

Page 1



### Product description

Order number	47.409.001
Product description	Sample tube, HbA1c determination from capillary blood, nominal volume: 1,000 µl, cap black, preparation: Roche HbA1c haemolysis solution, screw cap, (LxØ): 82 x 11.5 mm, flat base, with paper label, 50 piece(s)/bag

### Product characteristics

Type of collection	capillary
Type of preparation	Roche HbA1c haemolysis solution
Label/ Print	with paper label
Colour of print/label	black
Closure type	screw cap
Application	capillary blood collection
Example of use	HbA1c determination from capillary blood
Base shape	flat base
Volume of work	10 µl
Shape of inner tube	conical

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

Sample tube, Roche HbA1c haemolysis solution,  
1,000 µl, cap black, (LxØ): 82 x 11.5 mm, with  
paper label

Page 2

### Size

Sample volume	1.000.00 µl
Diameter	11.5 mm
Length of product	82 mm
Length including cap	82 mm
Length excluding cap	66 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	black
Tube material	Polypropylene (PP)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	37500
Depth of box	119 mm
Width of box	99 mm
Height of box	98 mm
Depth of case	513 mm
Width of case	293 mm
Height of case	140 mm
Case volume	0.021 cbm
Weight of product	0.0062 kg
Weight of case	3.71 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Sample tube, Roche HbA1c haemolysis solution,  
1,000 µl, cap black, (LxØ): 82 x 11.5 mm, with  
paper label**

**Page 3**

EAN of inner box                      4038917106382  
EAN case                                    4038917078399

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**