



Product description

Order number 95.963.007
Product description Wall holder, for Multi-Safe twin plus, material: stainless steel

Product characteristics

Example of use disposal of pointed and sharp objects
Suitable for Multi-Safe twin plus

Size

Width of product 186 mm
Length of product 165 mm

Material & colours

Product material stainless steel

Purity & certification

Product category no medical device | IVD
Latex-free yes
Batched no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	20
Piece(s) / pallet	2880
Depth of case	323 mm
Width of case	213 mm
Height of case	145 mm
Case volume	0.01 cbm
Weight of product	0.246 kg
Weight of case	5.24 kg
EAN of inner box	4038917900997
EAN case	4038917901000

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

