




**Product Specification**  
**Urine container, 90 ml, (ØxH): 60 x 65 mm, PP,  
with paper label, transparent**



### Product description

Order number	75.1435.004
Product description	Urine container, max. working volume: 90 ml, (ØxH): 60 x 65 mm, Ø opening: 65 mm, transparent, graduated, material: PP, screw cap, with paper label, cap assembled, sterile, 400 piece(s)/case

### Product characteristics

Label/ Print	with paper label
Graduation	yes
Closure type	screw cap
Volume of work	90 ml

### Size

Sample volume	90 ml
Max. volume	90 ml
Diameter	60 mm
Diameter of opening	65 mm
Height of product	65 mm
Length excluding cap	60 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine container, 90 ml, (ØxH): 60 x 65 mm, PP,  
with paper label, transparent**

**Page 2**

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

**Packaging**

Minimum order qty. 400  
Type of smallest case  
subpackaging  
Piece(s) / inner box 400  
Piece (s) / outer case 400  
Piece(s) / pallet 7200  
Case volume 0.1062 cbm  
Weight of product 0.0201 kg  
Weight of case 8.23 kg  
EAN of inner box 4038917385220  
EAN case 4038917385213

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**