



Product description

| | |
|---------------------|--|
| Order number | 75.1437.100 |
| Product description | Urine container, max. working volume: 75 ml, (ØxH): 65 x 48 mm, transparent, graduated, material: PP |

Product characteristics

| | |
|----------------|-------|
| Graduation | yes |
| Volume of work | 75 ml |

Size

| | |
|-------------------|-------|
| Sample volume | 75 ml |
| Max. volume | 75 ml |
| Diameter | 65 mm |
| Height of product | 48 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Urine container, 75 ml, (ØxH): 65 x 48 mm, PP,
transparent**

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1000
Piece(s) / inner box 25
Piece (s) / outer case 1000
Piece(s) / pallet 28000
Case volume 0.065 cbm
Weight of product 0.0092 kg
Weight of case 9.65 kg
EAN of inner box 4038917129138
EAN case 4038917074711

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT