



Product Specification

Multi-purpose container, max. volume: 250 ml,
(LxØ): 100 x 73 mm, PS, white

Page 1



Product description

Order number	75.560
Product description	Multi-purpose container, max. volume: 250 ml, (LxØ): 100 x 73 mm, Ø opening: 73 mm, white, material: PS, 50 piece(s)/bag

Product characteristics

Graduation	no
Application	collection and storage of urine
Volume of work	250 ml

Size

Sample volume	250 ml
Max. volume	250 ml
Diameter	73 mm
Diameter of opening	73 mm
Height of product	100 mm
Length of product	100 mm
Length excluding cap	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-22

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Multi-purpose container, max. volume: 250 ml,
(LxØ): 100 x 73 mm, PS, white

Material & colours

Product material Polystyrene (PS)
Colour of product white

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 500
Type of smallest
subpackaging bag
Piece(s) / inner box 500
Piece (s) / outer case 500
Piece(s) / pallet 20000
Depth of case 598 mm
Width of case 398 mm
Height of case 165 mm
Case volume 0.0393 cbm
Weight of product 0.0076 kg
Weight of case 4.3 kg
EAN of inner box 4038917339902
EAN case 4038917339896

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

