



## Product Specification

Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, PP, with safety label, with integrated transfer unit, transparent

Page 1



### Product description

Order number	75.562.400
Product description	Container with screw cap, collection and storage of urine, max. working volume: 100 ml, (ØxH): 57 x 76 mm, Ø opening: 62 mm, transparent, graduated, material: PP, screw cap, cap: yellow, with safety label, with integrated transfer unit, cap assembled, sterile, 5 piece(s)/tube bag

### Product characteristics

Label/ Print	with safety label
Colour of print/label	white/red/black
Graduation	yes
Closure type	screw cap
Design	with integrated transfer unit
Application	collection and storage of urine
Example of use	closed urine transfer into a V-Monovette® Urine
Volume of work	20 ml - 100 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, PP, with safety label, with integrated transfer unit, transparent

Page 2

### Size

Sample volume	100 ml
Max. volume	100 ml
Diameter	57 mm
Diameter of opening	62 mm
Height of product	76 mm
Length excluding cap	72 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow

### Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### Packaging

Minimum order qty.	200
Type of smallest subpackaging	tube bag
Piece(s) / inner box	200
Piece (s) / outer case	200
Piece(s) / pallet	5600
Depth of case	588 mm
Width of case	398 mm
Height of case	255 mm
Case volume	0.0597 cbm
Weight of product	0.0173 kg
Weight of case	4.31 kg
EAN of inner box	4038917236805
EAN case	4038917236799

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**