



## Product Specification

Multi-purpose container, max. volume: 125 ml,  
(LxØ): 85 x 62 mm, graduated, PP, transparent

Page 1



### Product description

Order number	75.565
Product description	Multi-purpose container, max. volume: 125 ml, (LxØ): 85 x 62 mm, Ø opening: 62 mm, transparent, graduated, material: PP, 50 piece(s)/bag

### Product characteristics

Graduation	yes
Application	collection and storage of urine
Volume of work	125 ml

### Size

Sample volume	125 ml
Max. volume	125 ml
Diameter	62 mm
Diameter of opening	62 mm
Height of product	85 mm
Length of product	85 mm
Length excluding cap	85 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-22

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Multi-purpose container, max. volume: 125 ml,  
(LxØ): 85 x 62 mm, graduated, PP, transparent**

**Page 2**

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

### Purity & certification

Product category	no medical device   IVD
Batched	yes

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	18000
Depth of case	578 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0415 cbm
Weight of product	0.0107 kg
Weight of case	5.93 kg
EAN of inner box	4038917082525
EAN case	4038917001601

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-22

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**